



BioAnalyt

Rapid Measurement of **Iron** and **Vitamin A** in Flour with iCheck



Special Report, September 2020
Revised on January 2023



Fortification of Wheat and Maize Flours

Wheat, maize, and rice are the most widely consumed grains around the world, present in about 95% of all diets. However, modern processing of these grains removes most of the bran and germ, the parts of the grains that are naturally rich in vitamins and minerals. People who rely on industrially processed wheat and maize flours or polished rice to supply most of their daily energy needs are, therefore, at a greater risk of developing micronutrient deficiencies.¹

Large-Scale Food Fortification (LSFF) is the addition of vitamins and minerals (e.g., iron, zinc, folic acid, vitamin A, and vitamin D) to staple foods such as wheat or maize flour, edible oils and fats, salt, sugar, and rice. LSFF is considered to be one of the most cost-effective public health interventions to increase the consumption of critical vitamins and minerals. Today cereal grains and flours are consumed at an average quantity of 200 grams/person/day, making these a suitable vehicle in many countries to deliver important micronutrients at the population level.²


In the early 1900s, deficiencies in B vitamins were prevalent in the USA as a result of the widespread consumption of industrially processed wheat flour. To tackle these deficiencies, mandatory fortification of processed flours was introduced in 1939 in the USA. At that time, wheat flour and white bread were enriched with thiamine (B1), riboflavin (B2), niacin (B3), iron, and calcium.³ A year later, in 1940, the United Kingdom also mandated fortification of wheat flour with thiamine. Denmark followed suit in 1953 by making the addition of thiamine, riboflavin, and iron to wheat flour mandatory.¹ Today, about 67 countries mandate wheat flour fortification, and 17 mandate both wheat and maize flour fortification.^{4,5} Both flours are commonly fortified with iron, zinc, B vitamins, and sometimes with vitamin A.

Vitamin and mineral premix for flour fortification costs an average of US \$4 per metric ton of flour for large-scale producers, which equals around US \$0.05-0.07 per person per year. The process of mixing in additives during flour production is well known to flour millers, since adding flour improvers is a common practice. Consequently, the addition of micronutrients to the production process carries a negligible cost once the equipment for fortification is put in place at a flour mill.^{1,5}

Furthermore, the benefits of fortifying flour with micronutrients far outweighs the cost of fortification. For every dollar spent on flour fortification, US \$84 is saved through reduction of anemia, a condition affecting one in three people globally.⁵ Sustained flour fortification has also been associated

with a 2.4% reduction in anemia prevalence among non-pregnant women of childbearing age, and an increase in serum retinol in populations from 0.94 to 1.06 $\mu\text{mol/L}$.^{6,7}

Fortification of flour has a noticeable effect on reducing neural tube defects in newborns, which are caused by insufficient folic acid stores. After Chile mandated flour fortification with folic acid in 2000, hospitals recorded a steady decrease in neural tube defects from 18.6 per ten thousand births in 1999 to 8.39 per ten thousand in 2009.⁵



LEGEND

| Mandatory cereal grain fortification legislation | |
|--|--|
| Grey | No mandatory legislation or data not available |
| Yellow | Wheat flour |
| Orange | Rice |
| Green | Wheat flour and maize flour |
| Light blue | Wheat flour and rice |
| Dark blue | Wheat flour, maize flour, and rice |

References:

- ¹ Hoogendoorn et al. 2FAS EC. 2016.
- ² Klemm et al. Food Nutr. Bull. 2010.
- ³ Bashai & Nalubola. EDCC. 2013

⁴ GFDx. 2020.

⁵ FFI, 2022.

⁶ Barkley et al. Br. J. Nutr 2015.

⁷ Dary, USAID. 2016.



How it works

Measuring Iron and Vitamin A in Flours with iCheck


To facilitate the monitoring and quality control of flour fortification, BioAnalyt has developed iCheck devices which enable the reliable measurement of the concentrations of vitamin A and iron in wheat and maize flours.

iCheck Iron is a portable, single-wavelength photometer that quantitatively measures iron in multiple food matrixes. Similar to reference methods, iCheck iron measures total iron in food samples (i.e., intrinsic iron as well as iron added as either ferrous fumarate, ferrous sulfate, NaFeEDTA, or ferric pyrophosphate).

Intrinsic iron is natural iron present in organic samples. In flour the intrinsic iron content may be between 10 and 60 mg/kg.


iCheck Fluoro is a portable, single-wavelength fluorometer that quantitatively measures vitamin A in foods and biological substances. iCheck Fluoro measures added vitamin A as retinyl palmitate and retinyl acetate.

IRON MEASUREMENT



* For ferrous fumarate, ferrous sulphate, and ferric pyrophosphate iron fortificants, it is recommended that your flour sample be diluted in 0.2M-0.4M HCl solution. These iron fortificants are partially soluble or insoluble in water. NaFeEDTA fortificant is soluble in water, hence water can be used as a diluent. It is recommended also to measure intrinsic iron in flour samples using 0.2M - 0.4M HCl

VITAMIN A MEASUREMENT



** When measuring vitamin A in flour, it is recommended to first measure an unfortified sample of the same flour to assess if the food matrix has innate fluorescence. This may cause iCheck Fluoro to display overestimated results.

Contact us to get **detailed protocol and training:**

Email: support@bioanalyt.com
WhatsApp: +49 3328 35150034



Performance of iCheck Iron and iCheck Fluoro are Comparable to Reference Methods

To ensure the reliability and accuracy of results, measurements obtained with iCheck devices are compared to results from traditional laboratory methodologies, such as High-Performance Liquid Chromatography (HPLC) for vitamin A and mass spectrometry (ICP-MS) for iron.

Wheat and maize flours were fortified with vitamin A at 1.0 - 10.0 mg retinol equivalents per kg (mg RE/kg) and measured with iCheck Fluoro and HPLC in two different studies. A good linear relationship was obtained between results from the two methods, and a correlation of $R^2=0.91$ to $R^2 = 0.97$ was observed.^{8,9} In another study, iCheck Fluoro results were compared to expected concentrations in wheat and maize flour samples spiked with vitamin A premix at 1 to 6 ppm. Recovery for wheat flour was at $116 \pm 2\%$, and for maize flour $113 \pm 5\%$.¹⁰

Results with iCheck Iron were compared to results with ICP-MS in wheat flour spiked with ferrous fumarate. The samples were diluted in 0.2M HCl for analysis with iCheck Iron, and the results showed excellent correlation with ICP-MS results ($R^2=0,96$), and a maximum coefficient of variation (CV) of 12% was observed. iCheck Iron results were compared to expected concentrations in wheat and maize flour samples spiked with NaFeEDTA at 15 to 90 ppm. The recovery of added iron for both types of flours was at $98 \pm 3\%$. A maximum CV of 6% was observed for wheat flour and 5% for maize flour, respectively. Recovery and precision was further assessed in wheat flour samples spiked with premix containing ferrous sulfate at 20 to 50 ppm iron. The recovery of iron was $101 \pm 7\%$, and maximum CV observed was 8%.


References:

⁸ Laillou et al, Food Nutr. Bull. 2014

⁹ Jaeger A, FST Mag. 2014

¹⁰ BioAnalyt internal data

Benefits of iCheck



- **Speed:** Results in 5 to 60 minutes.
- **Economy:** Cost is only 10% of conventional lab methods.
- **Easy implementation:** Only 1 day of training is required.
- **Scalability:** Portable, with no set-up calibration required.
- **Accuracy:** Performance is comparable to reference lab methods.

iChecks are manufactured in Germany, used in over 80 countries and validated against standard laboratory methods. Learn more at www.bioanalyt.com/products.

