

Quick Start Guide for iCheck Vitamin E

Whole Blood and Serum Samples

DOs



- Always wear gloves.
- Collect blood into tubes that contain anticoagulants.
- Slowly draw exactly 500 μ L into syringe.



- Ensure there are no air bubbles inside the syringe.
- Slowly inject 500 μ L through the red septum of the reagent vial.



- Shake the reagent vial with the sample vigorously for 10 seconds.
- Let the vial stand still for min 5 minutes.



- Two distinct phases must be observed inside the vial before measurement.
- Upper clear phase of 1 ± 0.2 cm high must be observed.



- Measure the vial with iCheck Vitamin E following instructions in the User Manual.

DON'Ts

- *Do not freeze whole blood samples.*

- *Do not recap needles to avoid injury.*
- *Do not dispose of used needles & syringes into common waste. Use specialized biohazard waste containers.*

- *Do not expose the vial with the sample to direct sunlight.*

- *Do not measure the vial if there are particles visible in the upper clear phase.*

- *Do not re-use the vial*
- *Do not dispose of used vials into common waste. Treat them as hazardous waste.*

Would you like to receive more information? Please contact us at

contact@bioanalyt.com • Tel.: +49 33 28 35 15 000 • WhatsApp: +49 332 835 150 034
BioAnalyt GmbH • Rheinstr. 17, 14513 Teltow, Germany • www.bioanalyt.com

